

UltraClean-DP-L Tabletop Ultrasonic Cleaning Machine

Features

- Automatic sweep function
- Degas function
- Ultrasonic power adjustment
- LCD control panel
- 1-99 minutes adjustable working time
- 1-80°C adjustable operating temperature
- Liquid drain system
- Complete stainless steel tanks and outer body
- High efficiency ultrasonic transducer system
- Stainless steel lid
- Stainless steel parts basket



LCD display



Degassing



Power adjustable



Heated



Timer adjustable

info@sterilmed-it.com

www.sterilmed-it.com

Address: Via Meravigli , Milano, Italy

Phone : +39 345 678 9012

UltraClean28-DP-L	
General specifications	
Tank dimensions (w*I*h)	505*300*200mm
Basket dimensions (w*I*h)	460*260*160mm
External dimensions (w*I*h)	520*320*400mm
Tank volume	28lt
Tank material	AISI 304 tank
Tank drain	¼" ball valve
Liquid level sensor	-
Heat isolation	Glass wool
Operating voltage	230V AC/50Hz-60Hz
Ultrasonic	
Ultrasonic power	600W
Ultrasonic transducers	10 pcs
Ultrasonic surface position	Bottom
Ultrasonic frequency	40 kHz
Heating power	600W
Control panel	LCD screen control