



Washer & Disinfector

Sterilmed Medical Washer Disinfectors available to wash and disinfect heat-resistant and reusable materials like Surgical Instruments, Operating Room Containers and Laboratory Instruments.

General Features

High Quality Stainless Steel (316 or 304 AISI)
Used in Hospitals, Outpatient Clinics, Rehabilitation Centers and Laboratories ,
Ergonomic design for easier and safer loading unloading operations,
3 water inlets for cold, hot and distilled water,
Easy one side access for service and maintenance,
Expanded capacity for detergent and neutralizer,
Easy to remove and clean internal filter.

Safe and Sound

With security measures that protect the operator and the device,
Password protection,
Emergency stop button,
Tempered Glass door resistant to impact and high temperatures,
Doors locks once the programs start running,
Unable to open both doors at once in Septic-Aseptic models.

Control and Observe

Every step on a handy 7" Touchscreen and PLC,
Full-automatic programs, adjustable for specific calibrations ,
User-Friendly easy to use UI,
Numeric and graphic display ,

Integrated 40-character thermal printer,
Visual, written and audial warning systems,
Data storage of up 1000 cycles,
USB port for transferring saved data to other media,
Communications via RS 232 and RS485 Ethernet Module.



Versatility

Of our devices makes room for a variety of additional features and accessories such as:



Turbo Heater



Condenser Cooling



Thermal Printer



Rack Identification System



Drain Pump



For a Green Future

For a Green Future

As Sterilmed Medical, we want to leave a green future for next generations and we manufacture our products with that sensibility; a zero waste mentality. All materials used in manufacturing our devices are recyclable or eco-friendly and we constantly work on increasing the efficiency of devices for less water and energy consumption.

Washer & Disinfector Equipments



DIN Baskets

Compatibility with ISO and DIN standard storage shelving systems.



Anesthesia Rack

Specially designed inner trolley that is compatible with Anesthesia Medical Equipment / Accessories.




Inner Trolley/ Rack

Grid type loading rack for ISO standard baskets.

Flexibility in height adjustment of the levels maximize load capacity.

Compatibility with ISO standard storage shelves systems.



Slipper Rack

Transport Trolley

Easy load system for stackable baskets and containers.

Compatibility with ISO and DIN standard storage shelving systems.

Availability of height adjustable and fixed height transfer trolley.





Quality Without Compromise

Is our motto and we manufacture our devices assuring their quality in compliance with;

Device Classification:
Class IIb, acc. To EC MDD (Annex II)

Medical Device Directive:
E93/42 EEC as amended by directive 2007/47/EC

Low Voltage Directive:
EN 60601-2-040

EMC Directive:
EN 61326-1

Washer Disinfector Equipment Directive:
EN 15883-1

Sterilization – Washer Disinfector – Large sterilizers:
EN 285

Quality Management System Requirements:
EN - ISO 9001

Medical Devices – Quality management systems –
Requirements for regulatory purposes: ISO 13485

Environmental Management Systems –
Requirements with guidance for use: ISO 14001

Capacity and Dimensions

Model	Chamber (WxHxD) - mm			Litres	Overall (WxHxD) - mm		
WD-8	625	680	600	250	980	1930	820
WD-10	625	680	600	300	980	1930	820
WD-12	625	680	850	350	980	1930	1070
WD-15	625	680	850	400	980	1930	1070
WD-18	625	680	850	450	980	1930	1070

* Single / Double Automatic Vertical Sliding, Tempered Glass Door System options are available.

As the main manufacturer of STERILMED Products, besides Standard STERILMED Model products and designs, visual changes can be made on a project-specific basis.



www.sterilmed-it.com



info@sterilmed-it.com
www.sterilmed-it.com

Address: Via Meravigli , Milano, Italy

Phone : +39 345 678 9012



Washer & Disinfector and Equipments



Address: Via Meravigli , Milano, Italy

info@sterilmed-it.com
www.sterilmed-it.com